

3870EA Large Capacity Automatic Sterilizer



3870EA

| Item No. | Chamber Dimensions Ø x L in(mm) | Volume | Overall Dimensions DxWxH in(mm) | Tray Dimensions DxWxH in(mm) | No. of Trays | Standard Unwrapped Cycle Time | Voltage(V) Freq.(Hz)* | Power(W) Current(A) | Shipping Weight lbs(kg) |
|-------------------|------------------------------------|-----------|------------------------------------|------------------------------------------------------------------|-----------------|-------------------------------------|--------------------------|------------------------|---------------------------------|
| 3870EA with stand | 15"x30" (384x760) | 22gal/85L | 34.5"x26"x20.6" (875x665x525) | 26.6"x11.3"x1" (675x286x25) 26.6"x13.8"x1" (675x350x25) | 2 | Cold: 32 mins Hot: 23 mins | 230V 50/60 | 3000W 13.6A | 278 (126) Including Stand |

NOTE: cycle times are approximate and may vary depending on instrument loads. *Model 3870EA requires a voltage between 220V and 235V for proper operation. When applicable a buck/boost transformer (0.5KVA) may be required.

Features

- Closed door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- A long life electro polished chamber and door
- New brighter and easier to read panel
- Program parameters can be customized and stored
- Supplied with a 3 shelf stand (30"Wx36"Dx34"H)
- 2 year parts and labor warranty and 10 year chamber warranty

Control Panel



1. Program Selection
2. Detailed Digital Display
3. Programmable Keys
4. System Indicators


Optional Internal Printer



Documents: date, time, temperature, pressure.



Clean & Simple™ Ultrasonic Cleaners and Tablets




Instrument Pouch Rack



Steam Distiller



Chamber Brite Autoclave Cleaner



Tuttnauer USA Co., Ltd.

25 Power Drive, Hauppauge, NY 11788 Tel: +800 624 5836,
+631 737 4850, Fax: +631 737 0720
E-mail: info@tuttnauerUSA.com, www.tuttnauerUSA.com

